



# 3Y series balances

**INNOVATIVE SOLUTIONS — NEW POSSIBILITIES** 

### 3Y series balances

Innovative solutions designed for 3Y series offer completely new possibilities, not only within operation but also when it comes to weighing results presentation.

Featuring 5,7" colour touchscreen display, 3Y series balances increase customers satisfaction. Great sunlight legibility of LCD, combined with clear menu arrangement, provides maximum comfort and ease of operation, significantly improving efficiency.

High precision touch-panel enables access to the most important functions, whereas build-in IR sensors facilitate user touch-free balance operation.

Multilingual, intuitive menu allows balance operators to personalize their profiles. Semi-automatic levelling system controls if the balance is correctly levelled, whereas ambient conditions are monitored by a dedicated built-in system ensuring appropriate conditions for the workplace.

Accuracy of measurements is provided by an automatic adjustment function. Every single weighing operation is recorded in an Alibi memory.

#### Main window

- A Information on current operating mode and profile
- Information on logged user
- Oate, time, information on connection type, battery state etc.
- Levelling status icon
- Measurement results
- (i) Minimum sample weight icon
- 6 Ambient conditions monitoring icon
- Configurable field for additional information
- Quick access bar





# Accuracy with no compromises


## Ambient conditions monitoring

The 3Y series offers a unique solution for monitoring of ambient conditions at the balance workplace.

The monitoring is performed in real time by means of internal, optionally external, sensors (THB module).


Depending on a balance type and sensors used, the system monitors temperature, humidity, atmospheric pressure and air density.






The system automatically informs about detection of change or about exceeding specified limits. The information is presented on the main screen of a balance. Thanks to AutoCal system, the user is warned about a need for balance adjustment. Each readout from a module is saved to balance memory.


## Auto-levelling Semi-automatic levelling system



## Air buoyancy compensation


Air buoyancy compensation aims at limiting an error occurring during weighment process performed on samples of density distinctly different from density of mass standard used for balance adjustment. Sensors of ambient conditions module help to calculate current air density, as a result the balance constantly corrects measurement result of a weighed sample.





## Minimum sample weight According to USP regulations

It is possible for the 3Y series to inform whether the minimum sample weight value has been exceeded for particular thresholds or not (for value greater than specified, measurements are performed with expected uncertainty of measurement). Such approach allows to






Upon touching the minimum sample weight icon, the user activates a window containing information on minimal sample weight value and the earliest date for its determination.

## Auto-test

## Ensuring the best quality of measurements

Auto-test function provides the user with possibility of manual confirmation of the performed measurements quality (Autotest GLP). Additionally the user may benefit from Autotest Filter function, which helps one to select such balance settings (appropriate for given ambient conditions) that make it work with the best possible repeatability.




# **Ergonomics and comfort of operation**

## **User profiles**


#### **Authorization and balance settings personalization**

Unique user, settings and authorization profiles will satisfy the user with flexible balance settings personalization. USB interface provides options for copying and moving user profiles from one 3Y series balance to another.



## Information system Complex databases

Complex databases are an outstanding feature of RADWAG balances, there are 18 secure databases for 3Y series balances. Particular databases size is dynamically shared within 1 GB memory flash. Option of databases import and export enables ease of databases management as well as their copying and archiving.




6 — RADWAG 2014

### **Sensors**

#### For touch-free operation

It is possible to assign an action to a particular IR sensor in order to perform touch-free activation of a given function (wide range of functions available).



## **Auto-open**

#### Automatic opening of draft shield door

For a greater comfort some balances of the 3Y series have been equipped with an automatic device facilitating opening and closing a draft shield door.





## **Portability**

## Wireless communication between a balance and a terminal

Use of wireless communication allows to place the terminal up to 10 meters away from the weighing chamber. Battery-powered terminal can be operated for up to 8 hours. Cable connection of the weighing chamber and the terminal is an alternative option.




RADWAG 2014 — 7

# Data management and safety

## **Protecting data**User authorization levels

When there is a need for one balance to be operated by several users, the option of customizing access rights for particular functions may turn out to be indispensable. Four access levels ensure many possibilities of supervision over users and provide important data protection (e.g. formulas).




Numerous operations such as defining language of the menu, selecting a desired working mode or personalizing main screen layout may be limited according to the access level, with respective password protection.



Possibility of associating a given profile with a particular user allows such balance personalization, that upon log in, a given working mode and filters are automatically selected.

## Alibi Secure data storage partition




3Y series balances feature secure partition for data storing where all weighment data, reports, ambient conditions measurements are recorded and secured for a specified amount of time. All these can be easily restored in case there was such a need.

## Data archiving and exchange

The 3Y series offers complex archiving of databases, user profiles and data stored in Alibi memory. All the data can be exported, imported, copied and transferred between balances.





**8** — RADWAG 2014

# Reports and printouts

## Reports Reports database

Upon completion of each process, a respective report is generated and recorded in a proper report database.

The users have possibility to preview, print, export or archive reports but not only. They can also freely configure them.

Weighing
Date 2014.08.19
Time 14:48:50
Balance type XA 3Y
Balance ID 392543
Level status Yes
Product Calcium
Net 0.7502 g
Tare 24.23788 g
Gross 24.9881 g
Signature

## **Printouts**Flexibility of configuration

There are two printout types for 3Y series balances: standard (generated according to a fixed pattern) and non-standard, customized ones.

Standard printout comprises three sections: header [A], weighment data [B] and footer [C].


Each section can be freely adjusted by a user, it can also be extended with a nonstandard printout.


	Weighing	
	Date	2014.04.02
	Time	14:07:43
	Balance ID	419036
A	Operator	Admin
	Level status	Yes
	Product	Calcium
L	Packaging	Blister
Г	Temperature during mesaurements:	26.79 ℃
	Humidity during measurements:	24 %
	Pressure during measurements:	994 hPa
	Net	0.1118376 q
	Tare	0.1118378 g 0.5000000 q
	Gross	0.6118376 q
	Supplementary unit	0.5591880 ct
	Minimum sample status	0K
B	Net	0.1118071 g
	Tare	0.5000000 g
	Gross	0.6118071 g
	Supplementary unit	0.5590355 ct
	Minimum sample status	OK
	Net	0.1118071 g
	Tare	0.5000000 g
	Gross	0.6118071 g
	Supplementary unit	0.5590355 ct
L	Minimum sample status	OK
	Cinnakum	
Y	Signature	
L		

#### **Charts**

#### **Measurement visualization and statistics**

Selected working modes (Statistics, SQC), apart from generating report, offer possibility of creating a chart for a particular completed weighment. The balance allows to generate weighments charts (with mean value calculated) and probability distribution chart out of series of measurements. Each chart can be freely scaled, printed or saved to BMP file.





# Operations to be carried out on printouts

#### Import, export and printing

Technology used for 3Y series enables free exchange of printouts and labels (TXT or LB file format) between balances. The balances are compatible with vast range of PCL printers and label printers, wherein connection between printers and balances is established by RS 232, USB, Ethernet.

# Technical specification

	Ultra and microbalances	microbalances Analytical balances		
	UYA 3Y, MYA 3Y	XA 3Y	XA 3Y.A	XA 3Y.F
Verification M	-	•	•	•
Max capacity [Max]	0,8 - 21 g	52 – 510 g	52 – 310 g	52 g
Readability [d]	0,1 – 1 μg	0,01 – 0,1 mg	0,01 – 0,1 mg	0,01 mg
Minimum weight USP	0,8 - 3,2 mg	20 – 160 mg	20 – 160 mg	20 mg
Minimum weight ( $U = 1\%$ , $k = 2$ )	0,08 - 0,32 mg	2 – 16 mg	2 – 16 mg	2 mg
Pan Size	ø 16 mm ø 26 mm ø 100 mm ø 160 mm	ø 85 mm ø 100 mm	ø 85 mm ø 100 mm	ø 85 mm 210×254 mm
Stabilization time	5 s, (10 – 20 s for UYA 3Y)	3 – 5 s	3 – 5 s	5 s (30 s for filters)
Internal (automatic) calibration/adjustment	•	•	•	•
Display				5,7" colour,
Interfaces			2	2×USB, 2×RS 232,

### **Software**



#### **Database Editor**

Databases support for 3Y balances and PUE 7 and PUE HY terminals.



#### **PIPETTES**

Pipettes calibration and estimation of volume measurement errors in accordance with ISO 8655.



#### **E2R Weighing Records**

Record of weighments carried out by balances cooperating in a network.



#### **PW Win**

Balance and computer cooperation display of measurements and statistics.



#### **RAD Key**

Readout of balance data by means of defined Hot Key.



#### **E2R Moisture Analyzers**

Record of weighments carried out by moisture analyzers cooperating in a network.

	Precision balances						
AS 3Y	PS 3Y		APP 3Y	APP 3Y.1	MA 3Y		
	Town to	a sum					
•	•	•	•	-	-		
82 – 510 g	200 – 2000 g	1500 – 6000 g	10 – 35 kg	10 – 50 kg	60 – 200 g		
0,01 – 0,1 mg	1 – 10 mg	10 mg	0,01 - 0,1 g	0,01 - 0,1 g	0,1 – 1 mg		
20 mg	-	-	-	-	-		
2 mg	-	-	-	-	-		
ø 70 mm ø 85 mm	128×128 mm	195×195 mm	347×259 mm	347×259 mm	ø 90 mm		
3,5 – 6 s	1,5 – 2 s	1,5 s	3 – 4 s	3 – 4 s	-		
•	•	•	•	_	-		
resistive touchscreen display (640×480)							

Ethernet, 4 Inputs / 4 Outputs (digital), Wi-Fi \*

## **Optional accessories**

Anti-vibration tables, Thermal and dot matrix printers, Adapter for pipettes calibration, THB Ambient condition module, Automatic feeder, Barcode scanner, Laboratory ware holders, Density determination kit for liquids and solid bodies.

Complete offer is to be found on www.radwag.com website.

### **Services**

Calibration, Minimum sample weight estimation, Adjustment of weighing chamber temperature (moisture analyzers), PQ/OQ/IQ Validation, Periodic inspections.

11 RADWAG 2014

<sup>\*</sup> Wi-Fi interface is intended only for ultra-microbalances and microbalances.



