

Dewar Containers acc. to DIN EN ISO 16496, KGW-Isotherm

These Dewar containers are made of borosilicate glass 3.3 (DIN/ISO 3585). Extra tough design with blue plastic-coated cover made of galvanized iron or stainless steel. With thermally insulated plastic lid held by two springs. Detachable carry-handle. The lid has to have a 2 mm hole when the containers are used to transport and store liquid nitrogen or any other gassing refrigerants.

Cover made of blue plastic-coated and galvanized iron sheet

| Type | Vol. l | Ø A approx. mm | Ø B approx. mm | Ø C approx. mm | D approx. mm | E approx. mm | Order No. |
|------|-----------|----------------------|----------------------|----------------------|--------------------|--------------------|-------------|
| 26 | 1 | 100 | 120 | 132 | 150 | 215 | 17 50 30261 |
| 27 | 2 | 138 | 160 | 170 | 170 | 245 | 17 50 30271 |
| 28 | 3 | 138 | 160 | 170 | 230 | 305 | 17 50 30281 |
| 29 | 4 | 138 | 160 | 170 | 310 | 385 | 17 50 30291 |

Cover made of stainless steel

| Type | Vol. l | Ø A approx. mm | Ø B approx. mm | Ø C approx. mm | D approx. mm | E approx. mm | Order No. |
|------|-----------|----------------------|----------------------|----------------------|--------------------|--------------------|-------------|
| 26 | 1 | 100 | 120 | 132 | 150 | 215 | 17 50 30262 |
| 27 | 2 | 138 | 160 | 170 | 170 | 245 | 17 50 30272 |
| 28 | 3 | 138 | 160 | 170 | 230 | 305 | 17 50 30282 |
| 29 | 4 | 138 | 160 | 170 | 310 | 385 | 17 50 30292 |

