

ECO-740 LAB

Oil less vacuum pump



SPECIFICATIONS	
N° of cylinders:	2
Air flow	24 LPM
Max Vacuum	740 mmHg/ 29,1”Hg
RPM	1750
Motor Power	85W – Thermal protector
Net weight	4,15Kg
Maximum dimensions in mm:	LxWxH (255x196x198)
Voltage Range	220v/60Hz 220v/50Hz 127v/60Hz
IP-20	
High Chemical Resistant	Optional

- Performance in C.A 60Hz at sea level