

ECO-290 LAB

Oil less vacuum pump



SPECIFICATIONS	
N° of cylinders:	2
Air flow	55 LPM
Max Vacuum	740 mmHg/29,1”Hg
RPM	3450
Motor Power	100W- Thermal protector
Net weight	5,15Kg
Maximum dimensions in mm:	LxWxH (275x196x199)
Voltage Range	220v/60Hz 220v/50Hz 127v/60Hz
IP-20	
High Chemical Resistant	Optional

- Performance in C.A 60Hz at sea level