

**Photos:**

**Spindles**



No.1



No.2



No.3



No.4



No.0: test 1-100 mpa.s





Inner padding



Out carton packing



**Specifications:**

Model	RM-1 (NDJ-1)	RM-2 (NDJ-4)
Measurement Range	1-100,000 mpa.s	1-2,000,000 mpa.s
Spindles	No. 0, 1, 2, 3, 4	
Rotating Speeds	6 / 12 / 30 / 60 rpm/min.	0.3 / 0.6 / 1.5 / 3 / 6 / 12 / 30 / 60 rpm/min.
Interface	No	
Measurement Tolerance	± 5%	
Power	AC 220V ± 10% 50 Hz ± 10%	
Working Condition	Temperature 5-35 degree Relative humidity no more than 80%	
Attachments	<ul style="list-style-type: none"><li>● No. 0, 1,2,3,4 Spindles<ul style="list-style-type: none"><li>● Adapter</li></ul></li><li>● Protection Frame</li><li>● Foundation Plate</li><li>● Up-Down Frame</li><li>● English Manual</li></ul>	
Outline (mm)	280×290×450 mm	
Packing (mm)	480×320×420 mm	
Net/Gross Weight	6 kg / 10kg	
H.S. Code	90278099	
Min. Order	1 pc by Air or 10 pcs by Ship	
Delivery Date (1 pc)	3 days After Receiving Payment	