


## 1.4 JHX Series - Micro Melting point Apparatus

Pictures	Name	Model	Temperature range	Detection method	Storage capacity (support)	Heating rate	Data interface	Fan	Can be equipped with printer
	Micro melting point instrument	JHX-4	Room temperature--350°C	Microscopic visual inspection		(0.1-10) °C/Min	RS232	NO	NO
		JHX-4Plus	Room temperature-400°C	Microscopic visual inspection				YES	YES
		JHX-5	Room temperature--350°C	Microscopic video				NO	NO
		JHX-5Plus	Room temperature--400°C	Microscopic video				YES	YES