

ACID DIGESTION BOMB



- **Volume:** 100 mL to 750 mL
- **Pressure Range:** Up to 200 bar
- **Temperature Range:** Up to 250°C
- **Material:** SS-316 with PTFE Liner and cap
- Microwaveable digestion bombs

AMAR EQUIPMENT PVT. LTD.

Valson Compound, LBS Road, Bhandup (W), Mumbai - 400078, INDIA.

+91-22- 6225 5000 www.amarequip.com sales@amarequip.com / export@amarequip.com